

Certificate of Completion
COVID-19 Testing Practices and Procedures Training
Trecia Billi Martin

Date / Time	Student Score
April 1, 2021 6:31 pm	100

Covid-19 Testing Practices and Procedures Training

#	Question	Correct Answer	Student Answer	Result	Points Awarded
1	Who should be tested for COVID-19?	Any patients showing the symptoms of fever, cough, or shortness of breath.	Any patients showing the symptoms of fever, cough, or shortness of breath.	Correct	10
2	Which of the following personal protective equipment (PPE) should you use to prevent the spread of COVID-19?	Face shield, Gown, Disposable mask, Gloves	Face shield, Gown, Disposable mask, Gloves	Correct	10
3	Which of the following is a step in proper handwashing techniques?	Keep fingers pointed downward so water does not roll onto arms.	Keep fingers pointed downward so water does not roll onto arms.	Correct	10
4	Which of the following statements about donning PPE is true?	The top of a mask should be pinched snugly over the bridge of the nose.	The top of a mask should be pinched snugly over the bridge of the nose.	Correct	10
5	Select all of the items that are in a diagnostic test kit.	Label to add patient information to the specimen, Viral transport tube(s), Synthetic fiber swab(s)	Label to add patient information to the specimen, Viral transport tube(s), Synthetic fiber swab(s)	Correct	10
6	Where should the patient ID, unique specimen ID, specimen type, date and time, and collector initials be recorded?	Directly on the transport tube label	Directly on the transport tube label	Correct	10
7	Select all true statements about nasopharyngeal (NP) swabbing.	Contact the back of the nasal cavity with the swab, and rotate over the surface.,	Contact the back of the nasal cavity with the swab, and rotate over the surface.,	Correct	10

		Insert the sterile swab straight back into the patient's nostril., Position the patient's head facing upward and tilted back.	Insert the sterile swab straight back into the patient's nostril., Position the patient's head facing upward and tilted back.		
8	Select all of the areas that the oropharyngeal swab should contact.	Tonsillar area, Back of the throat	Tonsillar area, Back of the throat	Correct	10
9	Place the following steps in the order you would remove your PPE.	1, Gown, 2, Gloves, 3, Face shield, 4, Mask	1, Gown, 2, Gloves, 3, Face shield, 4, Mask	Correct	10
10	Which of the following statements about ideal protocol for storage and transport of the samples in transport tubes is true?	Specimens should be stored between 2°C and 8°C and shipped with a cold pack within 72 hours.	Specimens should be stored between 2°C and 8°C and shipped with a cold pack within 72 hours.	Correct	10